

- SOP No:** AHP 80
- SUBJECT:** Magnetic resonance imaging (MRI) in live mice
- POLICY:** This procedure must be performed by an experienced operator.
- PRECAUTIONS:** Gloves, long-sleeved gown, closed in shoes.
- EQUIPMENT:** 700 MHz microimaging MR system  
Surface coil  
Anaesthetizing chamber  
Isoflurane  
Respiratory sensor  
26g needle  
1ml syringe  
Gastrolyte  
Heating pad
- PROCEDURE:**
1. Restrain mouse
  2. Induce anaesthesia with isoflurane at a dose of 4-5% in a closed anaesthetizing chamber then maintain at 1-2% by nose cone for the duration of the scan.
  3. Set up mouse on tooth hook of surface coil.
  4. Use padding and tape to secure respiratory sensor pad and to restrict movement of anaesthetised animal.
  5. Place surface coil containing mouse into the MRI, ensure temperature = 30 degrees C.
  6. Monitor respiration throughout experiment, maintain at 30-40resp/minute.
  7. The mouse is in the MRI for a period no longer than 3 hours.
  8. At the end of the scan, remove the mouse and stop isoflurane.
  9. Administer 1ml of gastrolyte (warmed to room temperature) by sub-cutaneous injection to prevent dehydration.
  10. Return mouse to home cage placed on heating pad.
  11. Monitor mouse for one hour to ensure full recovery.

## **RECOMMENDATION**

**DATE ISSUED:**                   **04.12.2008**

**REVISED:**

## **REFERENCES**

**Natt O, Watanabe T, Boretius S, Radulovic J, Frahm J and Michaelis T. (2002) *High-resolution 3D MRI of mouse brain reveals small cerebral structures in vivo. Journal of Neuroscience Methods 120:203-209.***